

# MEASURING GEOMETRIES

## MS-BV 1-1000

Measuring system using  
spindles in stainless steel 316 L\*

Stainless steel discs measuring system for viscosity measurements in a 150 mL beaker. This system allows rapid and economical measurements. This system could be used by First RM and can also be used with the Rheomat RM 100, RM 200.

NAME	REF	DIM. (mm)	VOL (mL)	VISCOSITY RANGE
	BV 1-100 AXIS	117102	-	-
	BV centring device	117202	-	-
	BV Disc n°1	117001	Ø 45	120
	BV Disc n°10	117010	Ø 40	120
	BV Disc n°100	117100	Ø 20	120
	BV 1000 Axis	117101	Ø 4	120
	150-ml glass beaker	117150	Ø 50-52	150
	MS TI Tube	118001	Ø 50	150

\* This system allows viscosity to be measured quickly and economically.